

Please complete all sections below and return to LGC: axiopt@lgcgroup.com

Purchase order no.: (compulsory)

PRODUCT INFORMATION

*easi-tab*TM are microorganisms presented in an easy-to-use lyophilized format as either individual tablets or vials. The microorganisms present in the *easi-tab*TM can be quickly resuscitated when required and are presented with a Statement of Measurement, showing the average level and the range using 95% confidence limits. Further information about the use of *easi-tab*TM can be found in the [instruction sheet](#).

The approximate levels of microorganisms in each product are given for guidance but may vary from batch to batch, further information regarding batch levels can be found on the Statement of Measurement for each product.

Level code	Target level (cfu/unit)	Range (cfu/unit)
E6	5,000,000	1,000,000 - 9,999,999
E5	500,000	100,000 - 999,999
E4	50,000	10,000 - 99,999
E3	5000	1000 - 9999
E2	500	100 - 999
E1	50	<100

CONDITIONS OF SUPPLY

The shelf life of the products is 24 months from production date. Please check individual Statements of Measurement for information regarding shelf-life. Materials are shipped with ice packs and should be stored at <-15°C on arrival at the customer's laboratory. Please refer to the [instruction sheet](#) for further details regarding storage and expiry.

ACCREDITATION

LGC is a UKAS accredited ISO 17034 reference material producer (29821) of microbiological reference materials. New batches of products (within our schedule of accreditation) are actively being produced on a rolling basis to replace stock produced before accreditation was awarded.

APPLICATION FORM

PRODUCTS AVAILABLE

MICROORGANISMS

Please indicate the number of packs required in the grid:

Organism	Product code LGCMIC	Target level (cfu/unit)	Format	Test microorganisms*	Source	ATCC	WDCM	Pack of 10 units
<i>Alicyclobacillus acidoterrestris</i>	RM73-S1-E1**	50	Tablet	Thermophilic acidophilic bacteria	NCIMB 13137	49025	-	
<i>Alicyclobacillus acidoterrestris</i>	RM73-S1-E2	500	Tablet	Thermophilic acidophilic bacteria	NCIMB 13137	49025	-	
<i>Aspergillus brasiliensis</i>	RM80-S1-E1**	50	Tablet	Osmophilic mould	NCPF 2275	16404	00053	
<i>Aspergillus brasiliensis</i>	RM80-S1-E4	50,000	Tablet	Mould	NCPF 2275	16404	00053	
<i>Bacillus cereus</i>	RM05-S1-E4	50,000	Tablet	<i>Bacillus cereus</i> ; <i>Bacillus</i> species	NCTC 7464	10876	-	
<i>Bacillus spizizenii (subtilis)</i>	RM23-S1-E4	50,000	Tablet	Aerobic spore count	NCTC 10400	6633	00003	
<i>Bacillus spizizenii (subtilis)</i>	RM23-S1-E6**	5,000,000	Tablet	Aerobic spore count	NCTC 10400	6633	00003	
<i>Candida albicans</i>	RM79-S1-E4	50,000	Tablet	Yeast	NCPF 3179	10231	00054	
<i>Clostridium perfringens</i>	RM06-S1-E4	50,000	Tablet	<i>Clostridium perfringens</i> / <i>Clostridium</i> species	NCTC 8237	13124	00007	
<i>Clostridium perfringens</i>	RM06-S2-E1	50	Tablet	<i>Clostridium perfringens</i> / Membrane Filtration	NCTC 8797	-	-	
<i>Clostridium perfringens</i>	RM06-S1-E3	5000	Tablet	<i>Clostridium perfringens</i> / <i>Clostridium</i> species	NCTC 8237	13124	00007	
<i>Clostridium perfringens</i>	RM06-S1-E5	500,000	Tablet	<i>Clostridium perfringens</i> / <i>Clostridium</i> species	NCTC 8237	13124	00007	
<i>Cronobacter sakazakii</i>	RM68-S1-E1	50	Tablet	<i>Cronobacter</i> species	NCTC 11467	29544	02214	
<i>Enterococcus faecalis</i>	RM01-S1-E1	50	Tablet	Enterococci / Membrane Filtration	NCTC 775	19433	00009	
<i>Enterococcus faecalis</i>	RM01-S1-E2	500	Tablet	Enterococci / Membrane Filtration	NCTC 775	19433	00009	

*Microorganisms enumerated to obtain the reference value mentioned on the *easi-tab*™ Statement of Measurement.

**Please note that these reference materials are not currently within the scope of LGC's UKAS accreditation.

APPLICATION FORM

Organism	Product code LGC MIC	Target level (cfu/unit)	Format	Test microorganisms*	Source	ATCC	WDCM	Pack of 10 units
<i>Enterococcus faecalis</i>	RM01-S1-E3	5000	Tablet	Total viable count at 22°C and 37°C	NCTC 775	19433	00009	
<i>Enterococcus faecalis</i>	RM01-S1-E4	50,000	Tablet	Total viable count	NCTC 775	19433	00009	
<i>Escherichia coli</i>	RM03-S1-E1	50	Vial	<i>Escherichia coli</i> Membrane Filtration	NCTC 9001	11775	00090	
<i>Escherichia coli</i>	RM03-S1-E4	50,000	Vial	<i>Escherichia coli</i>	NCTC 9001	11775	00090	
<i>Escherichia coli</i>	RM03-S3-E2	500	Vial	<i>Escherichia coli</i> Membrane Filtration	NCTC 12923	8739	00012	
<i>Escherichia coli</i>	RM03-S3-E3	5000	Vial	<i>Escherichia coli</i>	NCTC 12923	8739	00012	
<i>Escherichia coli</i>	RM03-S3-E4	50,000	Vial	<i>Escherichia coli</i>	NCTC 12923	8739	00012	
<i>Escherichia coli</i> O157:H7	RM25-S1-E1	50	Tablet	<i>Escherichia coli</i> O157:H7 (non-toxigenic)	NCTC 12900	700728	00014	
<i>Lactococcus lactis</i>	RM07-S1-E4	50,000	Tablet	Lactic acid bacteria	NCTC 6681	19435	00016	
<i>Legionella bozemanii</i>	RM48-S1-E4	50,000	Vial	<i>Legionella</i> species	NCTC 11368	33217	-	
<i>Legionella bozemanii</i>	RM48-S1-E6	5,000,000	Vial	<i>Legionella</i> species	NCTC 11368	33217	-	
<i>Legionella pneumophila</i> sg1	RM47-S1-E4	50,000	Vial	<i>Legionella pneumophila</i>	NCTC 11192	33152	00107	
<i>Legionella pneumophila</i> sg1	RM47-S1-E6	5,000,000	Vial	<i>Legionella pneumophila</i>	NCTC 11192	33152	00107	
<i>Listeria monocytogenes</i>	RM09-S1-E1	50	Tablet	<i>Listeria monocytogenes</i> <i>Listeria</i> species	NCTC 11994	-	00019	
<i>Listeria monocytogenes</i>	RM09-S1-E2	500	Tablet	<i>Listeria monocytogenes</i> <i>Listeria</i> species	NCTC 11994	-	00019	
<i>Listeria monocytogenes</i>	RM09-S1-E4	50,000	Tablet	<i>Listeria monocytogenes</i> <i>Listeria</i> species	NCTC 11994	-	00019	
<i>Penicillium chrysogenum</i>	RM22-S1-E4	50,000	Tablet	Mould	NCPF 2802	-	-	

*Microorganisms enumerated to obtain the reference value mentioned on the **easi-tab™** Statement of Measurement.

**Please note that these reference materials are not currently within the scope of LGC's UKAS accreditation.

APPLICATION FORM

Organism	Product code LGC MIC	Target level (cfu/unit)	Format	Test microorganisms*	Source	ATCC	WDCM	Pack of 10 units
<i>Pseudomonas aeruginosa</i>	RM17-S1-E1	50	Vial	<i>Pseudomonas aeruginosa</i> / Membrane Filtration	NCTC 10662	25668	00114	
<i>Pseudomonas aeruginosa</i>	RM17-S1-E2**	500	Vial	<i>Pseudomonas aeruginosa</i> / <i>Pseudomonas</i> species	NCTC 10662	25668	00114	
<i>Pseudomonas aeruginosa</i>	RM17-S1-E4	50,000	Vial	<i>Pseudomonas aeruginosa</i> / <i>Pseudomonas</i> species	NCTC 10662	25668	00114	
<i>Raoultella planticola</i>	RM34-S1-E1	50	Vial	Coliforms / Membrane Filtration	NCTC 9528	-	-	
<i>Raoultella planticola</i>	RM34-S1-E2	500	Vial	Coliforms / Membrane Filtration	NCTC 9528	-	-	
<i>Saccharomyces cerevisiae</i>	RM20-S1-E4	50,000	Tablet	Yeast	NCTC 10716	9763	00058	
<i>Salmonella</i> ser. Abony	RM08-S1-E1	50	Tablet	<i>Salmonella</i> (Group B)	NCTC 6017	BAA-2162	00029	
<i>Salmonella</i> ser. Typhimurium	RM08-S2-E1	50	Tablet	<i>Salmonella</i> (Group B)	NCTC 12023	14028	00031	
<i>Salmonella</i> ser. Tranoroa	RM08-S3-E1	50	Vial	<i>Salmonella</i> (Group O:55)	NCTC 10252	700148	-	
<i>Salmonella</i> ser. Tranoroa	RM08-S3-E2**	500	Vial	<i>Salmonella</i> (Group O:55)	NCTC 10252	700148	-	
<i>Staphylococcus aureus</i>	RM04-S1-E1	50	Vial	<i>Staphylococcus aureus</i> / Membrane Filtration	NCTC 6571	9144	00035	
<i>Staphylococcus aureus</i>	RM04-S1-E4	50,000	Vial	<i>Staphylococcus aureus</i> / <i>Staphylococcus</i> species	NCTC 6571	9144	00035	

*Microorganisms enumerated to obtain the reference value mentioned on the **easi-tab™** Statement of Measurement.

**Please note that these reference materials are not currently within the scope of LGC's UKAS accreditation.

Please note that LGC use source strains under licence and other strain references are given for information only.

NCTC = National Collection of Type Cultures (UK)

NCPF = National Collection of Pathogenic fungi (UK)

NCIMB = National Collection of Industrial and Marine bacteria (UK)

ATCC = American Type Culture Collection (USA)

WDCM = World Data Centre for Microorganisms. WDCM numbers in bold refer to strains that are listed in EN ISO 11133 standard (Microbiology of food, animal feeding stuffs and water - Preparation, production, storage and performance testing of culture media)

APPLICATION FORM

STERILE MATRICES

Please indicate the number of packs required in the grid:

Sample type	Suitable for	Product code LGCMIC	Format	Pack of 1 x 50g
Sterile skimmed milk powder	Milk & dairy testing	RM52**	50g	
Sterile oatmeal	Food testing	RM54**	50g	
Sterile cocoa powder	Confectionery testing	RM56**	50g	
Sterile dried diced chicken	Meat testing	RM58**	50g	

**Please note that these reference materials are not currently within the scope of LGC's UKAS accreditation.

SPECIAL ORDERS

BULK ORDER FOR SINGLE PRODUCTION LOTS

Enables Quality Manager to compare results from each of their laboratories

For your convenience and efficiency, we can offer a pre-allocated batch service.

How it works:

- You pre-order a specific quantity of packs and request that every pack comes from a single batch/lot.
- This quantity is pre-allocated to you. Depending on stock levels a new batch may be manufactured, please note this can take up to 12 weeks.

These can be:

- Either sent immediately, or
- You call off a number of packs as and when you need them.

This invokes our storage and management service where we organise safe storage in the freezer in Bury. There is a nominal material management charge.

CUSTOM BESPOKE PRODUCT SERVICE

If you have a requirement for an organism or level which is not available in our standard range of products listed above

Please send an enquiry to axiopt@lgcgroup.com stating the organism needed, the desired level (expressed as cfu per vial/tablet), the required quantity for the first batch, any regularity of the supply required and the list of countries to be delivered to and confirmation that each of those sites have all special local licences that are required to import and receive that organism. Note some countries require a customs form to be completed.

We will ascertain if an appropriate culture can be secured, the costs for the bespoke production and an estimated time.

Please note there will be a minimum order quantity (of 40 packs of a single type).

Products made using this service are not within our UKAS scope of accreditation

APPLICATION FORM

YOUR DETAILS

Please complete your details below:

Send Test Materials to:	Send Invoices to:
Contact name:	Contact name:
Department:	Department:
Company:	Company:
Address:	Address:
Town/City:	Town/City:
County:	County:
Post/Zip Code:	Post/Zip Code:
Country:	Country:
Tel:	Tel:
E-mail:	E-mail:
	VAT no:

Existing customers

If you would like to change your marketing preferences, please visit our [preference centre](#) or contact axiopt@lgcgroup.com

New customers

LGC would like to contact you by email, phone or post about its products, services and related areas or research. LGC does not sell contact details to external parties. You may unsubscribe at any time.

- ☐ Yes, I am happy to be contacted about LGC's products, services and research
- ☐ No, I do not require marketing information about LGC's products and services at this time

By signing this application form, you are agreeing to comply [by LGC's Standard Terms and Conditions for the Supply of Laboratory Quality Products and Services](#).

Furthermore, for test samples containing microorganisms with a biosafety level (BSL) of 1 or 2, you are confirming that:

- You are complying with applicable laws and regulations
- You shall ensure that laboratory facilities and expertise of staff are adequate for the safe handling and disposal of these organisms.

Name

Date

GENERAL NOTES

Please note:

- Please contact LGC for prices.
- An additional fee to cover courier charges. Please contact LGC for up-to-date charges.
- VAT (charged at the prevailing rate) is payable by all participants based in the UK.
- All courier charges set by LGC do not include local import charges, taxes (etc.), which will be covered by the customer.
- Do not send remittance with this form, you will be invoiced after the despatch of the materials.
- Please note, for large orders there may be a 12-week lead time if samples are required from the same batch.